


Band Box Sizes!

Description	Size	Order Number
Med. box	3/4 x 5 1/8 x 7 3/4	DDMBB 34
Lg. box	3/4 x 10 3/4 x 7 3/4	DDLBB 34
Lg. w/ inst. compartment	3/4 x 10 3/4 x 7 3/4	DDLBBI 34





Place 3-4 bands in each slot. When banding a patient, work directly out of the Band Box! If the band doesn't fit your patient, replace it in the band box, and try on another. When finished simply close the lid and sterilize the entire band box! NO MORE dirty piles of bands to sort through at the end of the day!

Stainless Steel

Cassettes are constructed of stain and rust resistant stainless steel, making sterilization of the entire cassette a breeze!

Sterilizable

Compared to the virtually impossible to sterilize cardboard and plastic band boxes commonly used, our stainless steel band boxes are a necessity.

Compatible

With all the band companies today, Steri-Source offers over 40 different patterns to suit your individual needs. Just tell us which band company you use, we'll have the perfect band box for that company.

Added Bonus

All Steri-Source band boxes include a FREE pair of Slim Line tweezers.

Color Coding

As always, we want to help your office be as efficient as possible. All our cassettes, cargo boxes and band boxes can be color coded.

Easy Identification

Each compartment in our band boxes is laser engraved to your band companies sizes. We want your experience with our band boxes to be the best!



Band Box

Stainless steel orthodontic band box, featuring double hinge lid, dual locking latches and laser engraving

Description:

Orthodontic stainless steel band box designed to ensure the proper sterilization of the entire band box. Eliminates the wasted time spent sorting and sterilizing individual bands. With the Steri-Source band box, you work directly out of the band box, when finished with the procedure the entire band box is sterilized in your dry heat, autoclave or flash sterilizer.

Intended Use:

Place 2 or 3 of each size band in your band box, when banding a patient, work directly out of the box. If a band that has been fitted is not used, replace it in the band box, and continue to fit the proper size band. When finished with the band box, replace the tweezers, close the box and run the entire band box through your sterilization system. Remember to restock the bands that you have used, to ensure ultimate efficiency.

Product Specifications:

Size: Medium: 3/4" x 5 1/8" x 7 3/4" Large: 3/4" x 5 1/8" x 10 3/4"

Compatible with Ultra sonic cleaners, autoclaves and dry heat.

Reorder Numbers:

Medium band box
Large band box
DDLBB34
Large with instrument comp.
DDLBBI34

Accessory Products:

Other sterilization products offered from *Steri-Source inc.* include: COTTRELL Good *Vibrations* multi-purpose ultrasonic cleaner rapidly remove organic and inorganic soils. *Steri-Source inc.* silicone instrument and cassette identification bands, efficiency made easy! **SPS Medical** steam *plus* integrators, spore tests and indicator tape, immediate peace of mind!

Band Company Patterns available

"A" Company, American Ortho, Class One, Dentaurun, GAC, Lancer, Ormco, Rocky Mountain Orthodontics, TP, Unitek. If your band company is not listed here, contact us at 1-877-783-7435, we may now have your pattern available.

Band Box Labels

76

67

Teflon labels are available at no extra cost to identify the type of bands you are storing in the band box. See below:

1 1	2 2	_ 3	_4 4_	_ 5
54 45	5 5 5	6 6	6	6
6 6	7 7	7	7	7 7 7

*Labels shown are examples of Upper Bands. Labels also available for Lower Bands