



Steri-Source Inc.

Sterilization Systems for Today's Orthodontist

Steri-Source Sterilization Flow



1 Simply place your non-sterile cassette into the waiting transport rack, now you are ready for your next patient! Steri-Source gets you right back with your patients, rather than tied up in the sterilization center!

2 Once the transport rack is full, it is ready to be lowered into the recessed ultrasonic! Steri-Source offers you custom sterilization design, we suggest using a recessed ultrasonic to save time, space and eliminate mess.



Sterile-To-Go!

3 Your cassettes are now ready to be pre-wrapped! Wrap your cassettes at the pre-wrap station and then place them back into the transport rack for sterilization.

4 Slide the transport rack with cassettes loaded into your sterilizer.... when sterile, your cassettes are ready to be placed in "Sterile-To-Go" already set-up, wrapped and ready for your next patient! *No individual sorting, scrubbing or set-ups are necessary with Steri-Source!*



Place cassettes into transport rack.



Lower transport rack into ultra sonic.



Pre-wrap or pouch cassettes.



Sterilize transport rack and cassettes.